




Acuity

FDR+



WORKFLOW INNOVATION

RadmediX 2510 Lance Drive, Kettering, OH 45409 USA

Tel: 844.723.6334 | Fax: 513.278.0150 | info@radmedix.com

© 2019 RadmediX Inc. All rights reserved.



FDR+ ADVANCED SYSTEM

Acuity FDR+ x-ray system with a standard configuration of floor mounted tube stand, vertical bucky stand and bucky table. This system features an advanced, 10" user-friendly touchscreen that provides a seamless operation of your radiographic system. It displays graphical tube angulation, SID values, and tube to bucky auto sync. With a range of advanced safety workflow and positioning features, Acuity FDR+ is a high performance x-ray system offering unparalleled price-performance ratio.

10" Intuitive Touchscreen

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



Premium Wall Stand

Motorized and counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning.

Overhead handgrips to promote proper upright exams.



AEC Optional

FDR+ KEY HIGHLIGHTS



Auto Sync Tracking with Wall & Table Bucky's



Rotational Column



6 Way Elevating Float Top Table



Wireless Pedal



AEC Optional

Emergency Stop

EFFORTLESS POSITIONING

Motorized auto-tracking enables a wide array of patient positioning.



INTUITIVE 10" TOUCHSCREEN



Elevated Table Specifications

- Lowering Height: 21 in
- Patient Load: 550 lbs
- Tabletop dimensions: 32 x 86 1/2 in
- Integrated Emergency Stop

FDR STANDARD SYSTEM

The Acuity FDR is also offered in our standard configuration. Offering the same versatility of the Acuity FDR+, but with our standard non-touchscreen display and advanced features. The Acuity FDR Standard offers the same exceptional digital imaging quality but with standard features at a budget friendly price point.

4" Console Display

Preset source to image distance with Auto Sync feature.

Real time display of SID and tube angulation.

Suitable for any clinical application.



AEC Optional

Premium Wall Stand

Manual counterbalance adjustment controls with electromagnetic locks.

Convenient handgrips to optimize patient positioning

Overhead handgrips to promote proper upright exams.



FDR STANDARD FEATURES



Rotational Column



4 Way Float Non Elevating Table



Built In Foot Pedal



AEC Optional

